




Occlutech AFR cardiac implant

Compatibility chart & ordering information



The Procedure Pack
INCLUDES
The Accessory Medical Devices
Occlutech Pistol Pusher



D1: DIAMETER OF THE FENESTRATION
D2: DIAMETERS OF THE DISCS
H: HEIGHT OF THE WAIST

Procedure Pack Art. No.	Ø D1 [mm]	Ø H [mm]	Ø D2 [mm]	ODS* Size [F]	Occlutech Pistol Pusher (OPP) Color
65AFR08MP	8	5	21	12	●
65AFR08LP	8	10	21	12	●
65AFR10MP	10	5	23	14	●
65AFR10LP	10	10	23	14	●

Product availability is subject to local regulatory clearances. Information shown is subject to change. Contact your sales representative for information.
* Occlutech Delivery Set (ODS), sold separately.

AFR
ATRIAL FLOW REGULATOR



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com

