




Occlutech PFO Occluder

Compatibility chart & ordering information



The Procedure Pack
INCLUDES
The Accessory Medical Devices
Occlutech Pistol Pusher



Ø LA DISC: LEFT ATRIAL DISC DIAMETER
Ø RA DISC: RIGHT ATRIAL DISC DIAMETER

Procedure Pack Art. No.	PFO Size (D) [mm]	Ø LA Disc [mm]	Ø RA Disc [mm]	Occlutech Pistol Pusher* Art. No. / Color	ODSv1** Size [F]	ODS III** Size [F]
LA Single Layer						
18PFO25SP	8<D≤13	23	25	55PP165 ●	9	9
LA Double Layer						
19PFO18DP	D≤8	16	18	55PP125 ●	7	7
19PFO25DP	8<D≤13	23	25	55PP165 ●	9	9
19PFO30DP	13<D≤15	27	30	55PP165 ●	9	9
19PFO35DP	15<D≤18	31	35	55PP210 ●	11	11

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.

* Color defines Pusher color.

** Figulla Flex II PFO occluder is compatible with Occlutech Delivery Set, ODSv1 (Article no.: 51DSXXX, 45°) and Occlutech Delivery Set III, ODS III (Article no.: 98DSXXX, 45°).



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com

