



# Occlutech PFO Occluder

Compatibility chart & ordering information



The Procedure Pack  
**INCLUDES**  
The Accessory Medical Devices  
Flex-Pusher II



Ø LA DISC: LEFT ATRIAL DISC DIAMETER  
Ø RA DISC: RIGHT ATRIAL DISC DIAMETER

Procedure Pack Art. No.	PFO Size (D) [mm]	Ø LA Disc [mm]	Ø RA Disc [mm]	Flex-Pusher II* Art. No. / Color	ODSv1** Size [F]	ODS III** Size [F]
<b>LA Single Layer</b>						
18PFO25SF	8<D≤13	23	25	51FP120 ●	9	9
<b>LA Double Layer</b>						
19PFO18DF	D≤8	16	18	51FP100 ●	7	7
19PFO25DF	8<D≤13	23	25	51FP120 ●	9	9
19PFO30DF	13<D≤15	27	30	51FP120 ●	9	9
19PFO35DF	15<D≤18	31	35	51FP160 ●	11	11

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.

\* Color defines Pusher color.

\*\* Figulla Flex II PFO occluder is compatible with Occlutech Delivery Set, ODSv1 (Article no.: 51DSXXX, 45°) and Occlutech Delivery Set III, ODS III (Article no.: 98DSXXX, 45°).



Occlutech International AB  
Landskronavägen 2, SE-252 32 Helsingborg, Sweden  
order@occlutech.com, [www.occlutech.com](http://www.occlutech.com)

