



Occlutech PLD Square

Compatibility chart & ordering information



The Procedure Pack
INCLUDES

The Accessory Medical Devices
Occlutech Pistol Pusher



D1: PROXIMAL DISC LENGTH
D2: DISTAL DISC LENGTH
D3: WAIST DIAMETER (ONLY FOR SQUARE W*)
H: HEIGHT OF THE WAIST (3 MM FOR ALL SIZES)

Procedure Pack Art. No.	Occlutech Pistol Pusher			Introducing System ¹	Cook-Cook ² Loader	Cook-Terumo ² Sheath	Cook-Terumo ² Loader	Cook-Terumo ² Sheath	ODS III** 45° 1100 mm
	[mm]	[mm]	[mm]						
Square W*	D1	D2	D3						
60PLD04WP	7.5	9	4	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	6F
60PLD05WP	8.5	10	5	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	6F
60PLD06WP	9.5	11	6	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	6F
60PLD07WP	10.5	12	7	●	RCF-7F	KCFW-8F	RCF-7F	RSC04-7F	7F
Square T*	Def S	D1	D2						
62PLD03TP	3	7	8	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	6F
62PLD05TP	5	8.5	10	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	6F
62PLD07TP	7	11	12	●	RCF-7F	KCFW-8F	RCF-7F	RSC04-7F	7F

Def S = Defect Size

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.

* W and T represent the types of connections between the discs. W stands for waist and T stands for twist.

1. Ordering Reference Numbers:

Loader:	French Size	Brand	Order No.	Ref No.	Length [cm]
	6F	Cook Performer (RCF)	G07337	RCF - 6.0 - 38 - J	13 cm
	7F	Cook Performer (RCF)	G07338	RCF - 7.0 - 38 - J	13 cm
	9F	Cook Performer (RCF)	G07564	RCF - 9.0 - 38 - J	13 cm
	10F	Cook Performer (RCF)	G07417	RCF - 10.0 - 38 - J	13 cm
Delivery Sheath:	French Size	Brand	Order No.	Ref No.	Length [cm]
	7F	Cook Ansel Flexor (KCFW)	G49044	KCFW - 7.0 - 35 - 110 - RB - HFANL1 - HC	110 cm
	7F	Terumo Destination (RSC04)	RSC04		90 cm
	8F	Cook Ansel Flexor (KCFW)	G49046	KCFW - 8.0 - 35 - 110 - RB - HFANL1 - HC	110 cm
	8F	Terumo Destination (54-89001)	54-89001		90 cm
	9F	Cook Performer (RCFW)	G09190	RCFW - 9.0 - 38 - 75 - RB - MTS	75 cm
	9F	Cook Performer (RCFW)	G08694	RCFW - 9.0 - 38 - 85 - RB - MTS	85 cm
	10F	Cook Performer (RCFW)	G08709	RCFW - 10.0 - 38 - 75 - RB - MTS	75 cm
	10F	Cook Performer (RCFW)	G08903	RCFW - 10.0 - 38 - 85 - RB - MTS	85 cm

2. Please note that Occlutech PLDs should not be retracted more than 3 times from the long sheath once they are deployed at the implantation site.

** Occlutech PLD is compatible with Occlutech Delivery Set III, ODS III (Article no.: 98DLXXX, 45°).




Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com




Occlutech PLD Rectangular

Compatibility chart & ordering information



The Procedure Pack
INCLUDES

The Accessory Medical Devices
Occlutech Pistol Pusher



AXB: WAIST LENGTH TO WIDTH (ONLY FOR RECTANGULAR W*)
D1XD3: PROXIMAL DISC LENGTH TO WIDTH
D2XD4: DISTAL DISC LENGTH TO WIDTH
H: HEIGHT OF THE WAIST (3 MM FOR ALL SIZES)

Procedure Pack Art. No.	[mm]	[mm]	[mm]	Occlutech Pistol Pusher (OPP) Color	Introducing System ¹				ODS III** 45° 1100 mm
					Cook-Cook ² Loader	Sheath	Cook-Terumo ² Loader	Sheath	
Rectangular W*	AxB	D1xD3	D2xD4						
61PLD04WP	4x2	8x6	9x7	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	6F
61PLD06WP	6x3	10x7	11x8	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	6F
61PLD08WP	8x4	12x8	14x10	●	RCF-7F	KCFW-8F	RCF-7F	RSC04-7F	7F
61PLD10WP	10x4	14x8	16x10	●	RCF-7F	KCFW-8F	RCF-7F	54-89001-8F	8F
61PLD12WP	12x5	16x9	17x10	●	RCF-9F	RCFW-9F	N/A	N/A	9F
61PLD14WP	14x6	19x11	21x13	●	RCF-9F	RCFW-9F	N/A	N/A	9F
61PLD16WP	16x8	20x12	23x14	●	RCF-10F	RCFW-10F	N/A	N/A	10F
61PLD18WP	18x10	22x14	24x16	●	RCF-10F	RCFW-10F	N/A	N/A	10F
Rectangular T*	Def S	D1xD3	D2xD4						
63PLD05TP	5x3	9x6	10x8	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	6F
63PLD07TP	7x4	11x8	12x9	●	RCF-6F	KCFW-7F	RCF-7F	RSC04-7F	7F
63PLD10TP	10x4	14x8	16x10	●	RCF-7F	KCFW-8F	RCF-7F	54-89001-8F	8F
63PLD12TP	12x5	16x9	17x11	●	RCF-9F	RCFW-9F	N/A	N/A	9F

Def S = Defect Size

Note: For the 9F and 10F size Paravalvular Leak Devices there is no compatible Terumo Long Sheath.

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.

* W and T represent the types of connections between the discs. W stands for waist and T stands for twist.

1. Ordering Reference Numbers:

Loader:	French Size	Brand	Order No.	Ref No.	Length [cm]
	6F	Cook Performer (RCF)	G07337	RCF - 6.0 - 38 - J	13 cm
	7F	Cook Performer (RCF)	G07338	RCF - 7.0 - 38 - J	13 cm
	9F	Cook Performer (RCF)	G07564	RCF - 9.0 - 38 - J	13 cm
	10F	Cook Performer (RCF)	G07417	RCF - 10.0 - 38 - J	13 cm
Delivery Sheath:	French Size	Brand	Order No.	Ref No.	Length [cm]
	7F	Cook Ansel Flexor (KCFW)	G49044	KCFW - 7.0 - 35 - 110 - RB - HFANL1 - HC	110 cm
	7F	Terumo Destination (RSC04)	RSC04		90 cm
	8F	Cook Ansel Flexor (KCFW)	G49046	KCFW - 8.0 - 35 - 110 - RB - HFANL1 - HC	110 cm
	8F	Terumo Destination (54-89001)	54-89001		90 cm
	9F	Cook Performer (RCFW)	G09190	RCFW - 9.0 - 38 - 75 - RB - MTS	75 cm
	9F	Cook Performer (RCFW)	G08694	RCFW - 9.0 - 38 - 85 - RB - MTS	85 cm
	10F	Cook Performer (RCFW)	G08709	RCFW - 10.0 - 38 - 75 - RB - MTS	75 cm
	10F	Cook Performer (RCFW)	G08903	RCFW - 10.0 - 38 - 85 - RB - MTS	85 cm

2. Please note that Occlutech PLDs should not be retracted more than 3 times from the long sheath once they are deployed at the implantation site.

** Occlutech PLD is compatible with Occlutech Delivery Set III, ODS III (Article no.: 98DLXXX, 45°).



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com

