



Occlutech ASD Occluder

Compatibility chart

ASD




Ø LA DISC: LEFT ATRIAL DISC DIAMETER
 Ø WAIST: MAX. WAIST DIAMETER
 Ø RA DISC: RIGHT ATRIAL DISC DIAMETER

Figulla Flex II ASD Occluder Ref. No.	AS Defect Size (D) [mm]	Ø Waist [mm]	Ø LA Disc [mm]	Ø RA Disc [mm]	Pusher System Art. No. / Color		ODSv1* Size [F]	ODS III* Size / Color [F]
					Flex Pusher II	Pistol Pusher		
29ASD04	D≤4	4	11	9	51FP060 ●	55PP090 ●	7	6 ●
29ASD05	4<D≤5	5	14	11	51FP060 ●	55PP090 ●	7	6 ●
29ASD06	5<D≤6	6	16.5	12.5	51FP100 ●	55PP125 ●	7	6 ●
29ASD07	6<D≤7.5	7.5	18	14	51FP100 ●	55PP125 ●	7	6 ●
29ASD09	7.5<D≤9	9	20.5	16.5	51FP100 ●	55PP125 ●	7	6 ●
29ASD10	9<D≤10.5	10.5	22	18	51FP100 ●	55PP125 ●	7	6 ●
29ASD12	10.5<D≤12	12	27	23	51FP120 ●	55PP165 ●	9	8 ●
29ASD13	12<D≤13.5	13.5	28.5	24.5	51FP120 ●	55PP165 ●	9	8 ●
29ASD15	12<D≤15	15	30	26	51FP120 ●	55PP165 ●	9	8 ●
29ASD16	15<D≤16.5	16.5	31.5	27.5	51FP120 ●	55PP165 ●	9	8 ●
29ASD18	15<D≤18	18	33	29	51FP120 ●	55PP165 ●	9	8 ●
29ASD19	16.5<D≤19.5	19.5	34.5	30.5	51FP120 ●	55PP165 ●	10	9 ●
29ASD21	18<D≤21	21	36	32	51FP150 ●	55PP185 ●	11	10 ●
29ASD24	21<D≤24	24	39	35	51FP150 ●	55PP185 ●	11	10 ●
29ASD27	24<D≤27	27	42	38	51FP160 ●	55PP210 ●	12	11 ●
29ASD30	27<D≤30	30	45	41	51FP160 ●	55PP210 ●	12	11 ●
29ASD33	30<D≤33	33	48	43	51FP160 ●	55PP210 ●	12	12 ●
29ASD36	33<D≤36	36	52	46	51FP160 ●	55PP210 ●	12	14 ●
29ASD39	36<D≤39	39	54	49	51FP160 ●	55PP210 ●	12	14 ●
29ASD40	39<D≤40	40	55	50	51FP160 ●	55PP210 ●	12	14 ●

Individual product availability is subject to local regulatory clearance and may not be available in all territories.

* Figulla Flex II ASD occluder is compatible with Occlutech Delivery Set, ODSv1 (Article no.: 51DSXXX, 45*) and Occlutech Delivery Set III, ODS III (Article no.: 98DSXXX, 45*).



Occlutech International AB
 Landskronavägen 2, SE-252 32 Helsingborg, Sweden
 order@occlutech.com, www.occlutech.com

