




# Occlutech mVSD Occluder

Compatibility chart & ordering information




The Procedure Pack  
**INCLUDES**  
The Accessory Medical Devices  
Occlutech Pistol Pusher



D1: DIAMETER OF THE WAIST  
D2: DIAMETER OF THE DISCS  
H: HEIGHT OF THE WAIST (7 MM FOR ALL SIZES)

Procedure Pack Art. No.	Defect Size [mm]	Height [mm]	Ø D1 [mm]	Ø D2 [mm]	ODS* Size [F]	Occlutech Pistol Pusher (OPP) Color
71VSD04P	2<D≤4	7	4	10	6	●
71VSD06P	4<D≤6	7	6	13	6	●
71VSD08P	6<D≤8	7	8	15	7	●
71VSD10P	8<D≤10	7	10	17	8	●
71VSD12P	10<D≤12	7	12	19	8	●
71VSD14P	12<D≤14	7	14	22	8	●
71VSD16P	14<D≤16	7	16	24	9	●
71VSD18P	16<D≤18	7	18	26	10	●
71VSD20P	18<D≤20	7	20	28	11	●

Availability subject to regulatory approval. Product data may be subject to change. Contact your Occlutech representative for information.  
\* Occlutech Delivery Set: ODSv1 (Article no.: 51USXXX, 45°, 80cm) and ODS III (Article no.: 98DSXXX, 45°, 79 cm and 98USXXX, 180°, 79 cm).



Occlutech International AB  
Landskronavägen 2, SE-252 32 Helsingborg, Sweden  
order@occlutech.com, [www.occlutech.com](http://www.occlutech.com)

