


Product specifications & ordering information



D1: DISC DIAMETER
D2: WAIST DIAMETER
D3: SHANK DIAMETER
L: DEVICE LENGTH
AO: AORTA
PA: PULMONARY ARTERY

Occlutech PDA Occluder Ref. No.	L [mm]	ØD1 [mm]	ØD2 [mm]	ØD3 [mm]	ODSv1* 180°, 800 mm (Size)	ODS III* 180°, 790 mm (Size)	Occlutech Procedure Pack with Occlusions-Pusher Ref. No.
---------------------------------	--------	----------	----------	----------	----------------------------	------------------------------	--

STANDARD SHANK

42PDA05	4.25	9	3.5	5	51US006 (6F)	98US006 (6F)	42PDA05O
42PDA06	5.00	10	4	6	51US006 (6F)	98US006 (6F)	42PDA06O
42PDA07	6.05	11	5	7	51US006 (6F)	98US006 (6F)	42PDA07O
42PDA08	6.30	13	6	8	51US006 (6F)	98US006 (6F)	42PDA08O
42PDA10	7.00	16	8	10	51US007 (7F)	98US007 (7F)	42PDA10O
42PDA12	12.00	18	10	12	51US007 (7F)	98US007 (7F)	42PDA12O
42PDA15	14.00	20	12	15	51US008 (8F)	98US008 (8F)	42PDA15O
42PDA18	16.00	24	14	18	51US009 (9F)	98US009 (9F)	42PDA18O

LONG SHANK

43PDA05L	7.00	9	3.5	5	51US006 (6F)	98US006 (6F)	43PDA05LO
43PDA06L	7.50	10	4	6	51US006 (6F)	98US006 (6F)	43PDA06LO
43PDA07L	8.50	11	5	7	51US006 (6F)	98US006 (6F)	43PDA07LO
43PDA08L	9.00	13	6	8	51US006 (6F)	98US006 (6F)	43PDA08LO
43PDA10L	10.50	16	8	10	51US007 (7F)	98US007 (7F)	43PDA10LO

* Ordered separately.

Occlutech is a leading specialist provider of minimally invasive cardiac devices

Our mission is to improve the quality of life for people with heart conditions. The vision is to become a leading global specialist in cardiac devices, addressing congenital heart defects, stroke prevention and heart failure – because every beat counts for people living with heart conditions.

References:

- 1. Data on file
- 2. Steiner K, et al., Acta Paediatr. 2024 Apr;113(4):812-817
- 3. Godart F, et al., Int J Cardiol. 2018 Jun 15;261:54-57
- 4. Bilici M. et al., Pediatr Cardiol, 2017 Dec;38(8):1598-1605
- 5. Kudumula V. et al., J Invasive Cardiol. 2015 May;27(5):250-7



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com



Not approved or available for sale in the United States. Not available in all geographies. Please contact your Occlutech representative about the availability of Occlutech products in your area© Occlutech. All rights reserved. Occlutech is a registered trademark. MM_PDA_003_02




Occlutech PDA Occluder
Close with
confidence
a wide range of PDA
anatomies



Material and design



Unique and optimized braiding for high flexibility and adaptability



PET patches are sewn into the shank and in the disk for optimized closure and endothelialization promotion²




Occlutech Occlusions-Pusher



Close with confidence^{1,2,3}

Designed for closing a large variety of PDA morphologies²


- Wide range of diameter options (3.5 to 14 mm) and shank configurations to adapt to different anatomies²
- Flat aortic retention disc with no hub to minimize potential obstruction and metal protrusion in the aorta³
- High flexibility and adaptability with unique braiding structure²



Flat aortic retention disc


Designed for stability²

Thanks to its unique design with a wider pulmonary side and two shank length options²



Standard Shank

Long Shank




PDA defect


Occlutech PDA Occluder delivered

Safe, Effective, Proven^{2,3}

- High procedural success^{3,4}
- 100% closure rate at 18 months follow-up³
- Low adverse event rates^{2,3,5}



Occlutech PDA Occluder (42PDA06) in place.



Occlutech PDA Occluder (42PDA06) in place.