




Occlutech Sizing Balloon

# Defining the dimension




# Occlutech Sizing Balloon




## Defining the dimension

The Occlutech Sizing Balloon is intended to assist determining the size of a cardiac septal defect. Used in combination with radiographic or echocardiographic imaging, it helps in selecting the appropriate implant size for defect closure.



## Designations:

- (1) Usable length
- (2) X-ray markers
- (3) Soft tip
- (4) Inflation port
- (5) Catheter Shaft
- (6) Guide wire access



## Product specifications & ordering information

Article No.	Balloon Ø Nominal [mm]	Nominal Fill Vol. [ml]	Max. Fill Vol. [ml]	Balloon Ø at Max. Fill Vol. [mm]	Balloon Length [mm]	Catheter Shaft [F]	Usable Length [mm]	Compatible Guide Wire [inch]
53SB254	25	25	35	27	45	8	700	0.035
53SB355	35	60	90	40	50	8	700	0.035

## Occlutech is a leading specialist provider of minimally invasive cardiac devices

Our mission is to improve the quality of life for people with heart conditions. The vision is to become a leading global specialist in cardiac devices, addressing congenital heart defects, stroke prevention and heart failure – because every beat counts for people living with heart conditions.



**Occlutech International AB**  
Landskronavägen 2, SE-252 32 Helsingborg, Sweden  
order@occlutech.com, [www.occlutech.com](http://www.occlutech.com)

