




Occlutech Stiff Guide Wire

Guiding with confidence



Occlutech Stiff Guide Wire



Guiding with confidence

The Occlutech Stiff Guide Wire has been developed using a stainless-steel core and a PTFE coating. This combination gives well-balanced support for day-to-day cases, as well as advanced ones.

A soft, atraumatic J-tip, in combination with the torque control, allows navigation in anatomies with extreme dilations or severe angulation.¹

This stiff configuration coupled with the flexible tip allows a good push while minimizing risk of injuries.

- **Soft atraumatic J-tip** – supporting navigation
- **Stiff configuration** – supporting entry into the left atrium



Product specifications & ordering information

Article No.	Ø [inch]	Length [cm]	Tip shape [mm]
52GWR01	0.035	260	J - 3 mm

Reference:

1. *Instructions for Use – Guidewires for Central Circulatory Applications (EPflex), Rev. 1 (02/2024), REF 83026871*

Occlutech is a leading specialist provider of minimally invasive cardiac devices

Our mission is to improve the quality of life for people with heart conditions. The vision is to become a leading global specialist in cardiac devices, addressing congenital heart defects, stroke prevention and heart failure – because every beat counts for people living with heart conditions.



Occlutech International AB
Landskronavägen 2, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com



Not approved or available for sale in the United States. Product(s) may not be available in all markets. Please contact your Occlutech representative if you have questions about the availability of Occlutech products in your area. © Occlutech. All rights reserved. Occlutech is a registered trademark. MM_OSGW_001_01