




mVSD

Occlutech mVSD Occluder

Closing
muscular
VSDs



Product specifications & ordering information



D1: DIAMETER OF THE WAIST
D2: DIAMETER OF THE DISCS
H: HEIGHT OF THE WAIST (7 MM FOR ALL SIZES)

Occlutech mVSD Occluder Ref. No.	Defect Size [mm]	Ø D1 [mm]	Ø D2 [mm]	Height [mm]	ODSv1* 45° / 180°, 800 mm (Size)	ODS III* 45° / 180°, 790 mm (Size)	Occlutech Procedure Pack Flex Pusher II Ref. No.	Pistol Pusher Ref. No.
VSD04	2<D≤4	4	10	7	51US006	98DS006/98US006	71VSD04F	71VSD04P
VSD06	4<D≤6	6	13	7	51US006	98DS006/98US006	71VSD06F	71VSD06P
VSD08	6<D≤8	8	15	7	51DS007/51US007	98DS007/98US007	71VSD08F	71VSD08P
VSD10	8<D≤10	10	17	7	51DS008/51US008	98DS008/98US008	71VSD10F	71VSD10P
VSD12	10<D≤12	12	19	7	51DS008/51US008	98DS008/98US008	71VSD12F	71VSD12P
VSD14	12<D≤14	14	22	7	51DS008/51US008	98DS008/98US008	71VSD14F	71VSD14P
VSD16	14<D≤16	16	24	7	51DS009/51US009	98DS009/98US009	71VSD16F	71VSD16P
VSD18	16<D≤18	18	26	7	51DS010	98DS010/98US010	71VSD18F	71VSD18P
VSD20	18<D≤20	20	28	7	51DS011	98DS011/98US011	71VSD20F	71VSD20P

* Ordered separately.

Occlutech is a leading specialist provider of minimally invasive cardiac devices

Our mission is to improve the quality of life for people with heart conditions. The vision is to become a leading global specialist in cardiac devices, addressing congenital heart defects, stroke prevention and heart failure – because every beat counts for people living with heart conditions.



Occlutech International AB
Landskronavägen 9B, SE-252 32 Helsingborg, Sweden
order@occlutech.com, www.occlutech.com



Not approved or available for sale in the United States. Not available in all geographies. Please contact your Occlutech representative about the availability of Occlutech products in your area © Occlutech. All rights reserved. Occlutech is a registered trademark. MM_mVSD_001_02