
Occlutech PDA Occluder

Compatibility chart

PDA





D1:	DISC DIAMETER	L:	DEVICE LENGTH
D2:	WAIST DIAMETER	AO:	AORTA
D3:	SHANK DIAMETER	PA:	PULMONARY ARTERY

Occlutech PDA Occluder Ref. No.	**L** [mm]	**ØD1** [mm]	**Ø D2** [mm]	**Ø D3** [mm]	**ODSv1* 180°, 800 mm** (Size)	**ODS III* 180°, 790 mm** (Size)	**Occlutech PDA Occluder Procedure Pack Occlusions-Pusher** Ref. No.
STANDARD SHANK							
42PDA05	4.25	9	3.5	5	51US006 (6F)	98US006 (6F)	42PDA05O
42PDA06	5.00	10	4	6	51US006 (6F)	98US006 (6F)	42PDA06O
42PDA07	6.05	11	5	7	51US006 (6F)	98US006 (6F)	42PDA07O
42PDA08	6.30	13	6	8	51US006 (6F)	98US006 (6F)	42PDA08O
42PDA10	7.00	16	8	10	51US007 (7F)	98US007 (7F)	42PDA10O
42PDA12	12.00	18	10	12	51US007 (7F)	98US007 (7F)	42PDA12O
42PDA15	14.00	20	12	15	51US008 (8F)	98US008 (8F)	42PDA15O
42PDA18	16.00	24	14	18	51US009 (9F)	98US009 (9F)	42PDA18O
LONG SHANK							
43PDA05L	7.00	9	3.5	5	51US006 (6F)	98US006 (6F)	43PDA05LO
43PDA06L	7.50	10	4	6	51US006 (6F)	98US006 (6F)	43PDA06LO
43PDA07L	8.50	11	5	7	51US006 (6F)	98US006 (6F)	43PDA07LO
43PDA08L	9.00	13	6	8	51US006 (6F)	98US006 (6F)	43PDA08LO
43PDA10L	10.50	16	8	10	51US007 (7F)	98US007 (7F)	43PDA10LO

* Ordered separately



Occlutech International AB
Landskronavägen 9B, SE-252 32 Helsingborg, Sweden
order@occlutech.com, **www.occlutech.com**

Occlutech
Every beat counts

Not approved or available for sale in the United States. Not available in all geographies. Please contact your Occlutech representative about the availability of Occlutech products in your area. © Occlutech. All rights reserved. Occlutech is a registered trademark. MM_PDA_008_01