

Occlutech PFO Occluder

Compatibility chart

PFO





Ø LA DISC: LEFT ATRIAL DISC DIAMETER
Ø RA DISC: RIGHT ATRIAL DISC DIAMETER

Occlutech PFO Occluder	Defect Size (D)	Ø LA Disc	Ø RA Disc	ODSv1* 45°, 800 mm	ODS III* 45°, 790 mm	Occlutech PFO Occluder Procedure Pack	
						Flex Pusher II	Pistol Pusher
Ref. No.	[mm]	[mm]	[mm]	(Size)	(Size)	Ref. No.**	Ref. No.**
LA SINGLE LAYER							
18PFO25S	8<D≤13	23	25	51DS009 (9F)	98DS009 (9F)	18PFO25SF	18PFO25SP
LA DOUBLE LAYER							
19PFO18D	D≤8	16	18	51DS007 (7F)	98DS007 (7F)	19PFO18DF	19PFO18DP
19PFO25D	8<D≤13	23	25	51DS009 (9F)	98DS009 (9F)	19PFO25DF	19PFO25DP
19PFO30D	13<D≤15	27	30	51DS009 (9F)	98DS009 (9F)	19PFO30DF	19PFO30DP
19PFO35D	15<D≤18	31	35	51DS011 (11F)	98DS011 (11F)	19PFO35DF	19PFO35DP

* Occlutech PFO Occluder is compatible with Occlutech Delivery Set, ODSv1 (Article no.: 51DSXXX, 45°) and Occlutech Delivery Set III, ODS III (Article no.: 98DSXXX, 45°).
** Color defines Pusher color.



Occlutech International AB
Landskronavägen 9B, SE-252 32 Helsingborg, Sweden
order@occlutech.com, **www.occlutech.com**

Occlutech
Every beat counts

Not approved or available for sale in the United States. Product(s) may not be available in all markets. Please contact your Occlutech representative if you have questions about the availability of Occlutech products in your area. © Occlutech. All rights reserved. Occlutech is a registered trademark. MM_PFO_067_01