


Occlutech PLD Occluder

Compatibility chart

PLD








SQUARE:
D1: PROXIMAL DISC LENGTH
D2: DISTAL DISC LENGTH
D3: WAIST DIAMETER (ONLY FOR SQUARE W*)

RECTANGULAR:
AxB: WAIST LENGTH TO WIDTH (ONLY FOR RECTANGULAR W*)
D1xD3: PROXIMAL DISC LENGTH TO WIDTH
D2xD4: DISTAL DISC LENGTH TO WIDTH

GENERAL:
H: HEIGHT OF THE WAIST (3 mm FOR ALL SIZES)

Occlutech PLD Occluder[1]				**ODS III* 45°, 1100 mm**	**Occlutech PLD Occluder Procedure Pack****	
					Flex Pusher I	**Pistol Pusher**
Ref. No.	[mm]	[mm]	[mm]	(Size)	Ref. No.***	Ref. No.***
SQUARE WAIST	**D1**	**D2**	**D3**			
60PLD04W	7.5	9	4	98DL006 (6F)	60PLD04WF	60PLD04WP
60PLD05W	8.5	10	5	98DL006 (6F)	60PLD05WF	60PLD05WP
60PLD06W	9.5	11	6	98DL006 (6F)	60PLD06WF	60PLD06WP
60PLD07W	10.5	12	7	98DL007 (7F)	60PLD07WF	60PLD07WP
SQUARE TWIST	**D1**	**D2**	**Defect Size**			
62PLD03T	7	8	3	98DL006 (6F)	62PLD03TF	62PLD03TP
62PLD05T	8.5	10	5	98DL006 (6F)	62PLD05TF	62PLD05TP
62PLD07T	11	12	7	98DL007 (7F)	62PLD07TF	62PLD07TP
RECTANGULAR WAIST	**AxB**	**D1xD3**	**D2xD4**			
61PLD04W	4x2	8x6	9x7	98DL006 (6F)	61PLD04WF	61PLD04WP
61PLD06W	6x3	10x7	11x8	98DL006 (6F)	61PLD06WF	61PLD06WP
61PLD08W	8x4	12x8	14x10	98DL007 (7F)	61PLD08WF	61PLD08WP
61PLD10W	10x4	14x8	16x10	98DL008 (8F)	61PLD10WF	61PLD10WP
61PLD12W	12x5	16x9	17x10	98DL009 (9F)	61PLD12WF	61PLD12WP
61PLD14W	14x6	19x11	21x13	98DL009 (9F)	61PLD14WF	61PLD14WP
61PLD16W	16x8	20x12	23x14	98DL010 (10F)	61PLD16WF	61PLD16WP
61PLD18W	18x10	22x14	24x16	98DL010 (10F)	61PLD18WF	61PLD18WP
RECTANGULAR TWIST	**Defect Size**	**D1xD3**	**D2xD4**			
63PLD05T	5x3	9x6	10x8	98DL006 (6F)	63PLD05TF	63PLD05TP
63PLD07T	7x4	11x8	12x9	98DL007 (7F)	63PLD07TF	63PLD07TP
63PLD10T	10x4	14x8	16x10	98DL008 (8F)	63PLD10TF	63PLD10TP
63PLD12T	12x5	16x9	17x11	98DL009 (9F)	63PLD12TF	63PLD12TP

1. P17.05.002.07__Compatibility Chart PLD_ver07 * Ordered separately. ** Depending on availability in your country. *** Color defines Pusher color.



Occlutech International AB
Landskronavägen 9B, SE-252 32 Helsingborg, Sweden
order@occlutech.com, **www.occlutech.com**

Occlutech
Every beat counts

Not approved or available for sale in the United States. Not available in all geographies. Please contact your Occlutech representative about the availability of Occlutech products in your area. © Occlutech. All rights reserved. Occlutech is a registered trademark. MM_PLD_007_01